

# MUCILAIR™ VERSUS RAW 264.7 CELLS IN NANO- TOXICOLOGY

**TNO** triskelion bv

## OBJECTIVE

With increasing applications of engineered nano-materials (ENM) resulting in increasing exposure, the safety should be addressed along with the design of new materials. Animal testing may not be a desired test method for thousands of new nanomaterials because of the ethical concerns. Therefore, development of predictive *in vitro* toxicological screening can be valuable to rank ENM to determine priority for subsequent *in vivo* testing. In *in vitro* studies, nano-materials are predominantly studied in A549 or RAW 264.7 cell lines. However, using human 3D airway models opens up new possibilities in predictive *in vitro* testing of nanomaterials. These models consist of fully differentiated human epithelial cells and allow relevant exposure via air as they are cultured at an air-liquid interface. We compared the toxicity of SiO<sub>2</sub> and CeO<sub>2</sub> nanoparticles on MucilAir™ (EpiThelix Sarl) to RAW 264.7 macrophages.




FIGURE 1  
Schematic representation  
of the MucilAir™ model.

Frederique van Acker<sup>1</sup>, Mariska Gröllers<sup>2</sup>, Astrid Reus<sup>1</sup>, Yvonne Staal<sup>1</sup>, Esther Zondervan<sup>2</sup>, Frieke Kuper<sup>2</sup>


frederique.vanacker@tno.triskelion.nl

## METHODS

MucilAir™ inserts and RAW 264.7 cells were exposed for 24h to the nanoparticles via a droplet on the tissue surface and via the medium, respectively. Cytotoxicity was measured by LDH leakage (both) and TEER (MucilAir™) or MTT conversion (RAW 264.7). Various cytokines were analyzed in culture medium as a measure of inflammation. Oxidative stress and genotoxicity was evaluated by HO-1 expression and comet assay, respectively.

## RESULTS

In RAW 264.7 cells, SiO<sub>2</sub> and CeO<sub>2</sub> were cytotoxic at similar concentrations (figure 2b,c). A distinct induction of cytokine release (TNF-α only) was observed with SiO<sub>2</sub> only (figure 3b), whereas HO-1 expression and an increase in % tail DNA was only induced by CeO<sub>2</sub> (figures 4 and 5). In MucilAir™, no clear effects were observed for all endpoints with up to 10-fold higher concentrations (Figures 2a,b; 3-5).



FIGURES 2a and b  
LDH leakage and MTT data from RAW 264.7 cells treated for 24h with SiO<sub>2</sub> (fig 2a) and LDH leakage for MucilAir™ inserts treated for 24h with SiO<sub>2</sub>, CeO<sub>2</sub> (fig 2b).

Nano-sized SiO<sub>2</sub> and nano-sized CeO<sub>2</sub> give similar results in the LDH leakage and MTT assay (data CeO<sub>2</sub> not shown). IC<sub>50</sub> values are in the same order (data not shown). No cytotoxicity was observed in MucilAir™ based on TEER (data not shown) and LDH leakage for both SiO<sub>2</sub> and CeO<sub>2</sub>.




FIGURE 3a and b  
HO-1 expression data for RAW 264.7 cells (fig 3a) or MucilAir™ inserts (fig 3b) treated for 24h with SiO<sub>2</sub> or CeO<sub>2</sub>.

Although CeO<sub>2</sub> clearly induces HO-1 expression in RAW 264.7 cells, this is not observed in MucilAir™. SiO<sub>2</sub> does not induce HO-1 expression in either RAW 264.7 cells or MucilAir™.

1 TNO Triskelion, Utrechtseweg 48, Zeist, The Netherlands  
2 TNO, Utrechtseweg 48, Zeist, The Netherlands





FIGURE 4a and b  
TNF-α release for RAW 264.7 cells (fig 4a) or release of cytokines IL-8, IL-1α and TNF-α for MucilAir™ inserts treated for 24h with SiO<sub>2</sub> or CeO<sub>2</sub> (shown for RAW 264.7 cells only).

From a series of 13 cytokines investigated in RAW 264.7 cells, only TNF-α release was increased after exposure to SiO<sub>2</sub>. CeO<sub>2</sub> did not increase the release of TNF-α. In MucilAir™ neither SiO<sub>2</sub> nor CeO<sub>2</sub> demonstrated an increased release of IL-1α, IL-8 and TNF-α.



FIGURES 5a and b  
Comet assay for RAW 264.7 cells (fig 5a) or MucilAir™ inserts (fig 5b) treated for 24h with SiO<sub>2</sub> or CeO<sub>2</sub>.

CeO<sub>2</sub> induced a clear induction in the % tail DNA in RAW 264.7 cells. In MucilAir™ both damaged and undamaged cells were observed after exposure to CeO<sub>2</sub>, indicating that the nanoparticles might have an adverse effect on some cells. However, only a slight induction in % tail DNA was observed which needs further investigation. SiO<sub>2</sub> did not induce an increase in % tail DNA in either RAW 264.7 cells or MucilAir™.

## CONCLUSION

MucilAir™ appears less sensitive towards particle induced toxicity. As MucilAir™ enables a more relevant exposure, results from this assay could prove to be more realistic than experiments in cell lines. This will be investigated in further experiments.

## FUTURE PLANS

In future, we will optimize experimental conditions and further assess the applicability of 3D airway models, including MucilAir™, by different exposure systems and compare the results with both cell culture and *in vivo* inhalation data. Ultimately, these models may be useful in the first tier(s) of the safety evaluation of engineered nanomaterials.