


## Characterisation of IL-13 induced Airway Goblet Cell Metaplasia

Cindia Ferreira Lopes, Jimmy Vernaz, Nicolas Simonnet, Song Huang, Bernadett Boda, Ghislaine Arib, Samuel Constant



Epithelix Sàrl, 18 chemin des Aulx, CH-1228 Plan-les-Ouates, Geneva, Switzerland

### Introduction



Goblet cell metaplasia can be induced in a fully differentiated human airway epithelium following an IL-13 two-weeks treatment. Epithelix' previous data suggested that IL-13-induced epithelial remodeling was dependent on epithelium anatomical origin (Figure2) as bronchial tissue (MucilAir™-Bronchial) responded better to IL-13 than bronchiolar (SmallAir™) and nasal (MucilAir™-Nasal) tissues. To investigate this further, pooled versions of these epithelium models representing various anatomical regions were reconstituted using primary cell mix from several donors and were directly compared for their ability to acquire hallmark features of metaplasia following IL-13 exposure. Using both a newly developed platform based on a combination of immunofluorescence (IF) whole tissue staining and 3D confocal microscopy pictures analysis, and an IHC traditional approach, we demonstrated that MUC5AC increases in a dose-response manner in response to IL-13 concentrations ranging from 10 to 30 ng/ml in all 3 types of epithelia. Specifically, bronchial tissue gave the best response to IL-13 and nasal epithelium showed a less pronounced effect.



### Methods & Experimental Strategy



### Results



Measurements of TEER and CBF (Figure 5) indicated that epithelium barrier and ciliary functions remained unchanged compared to the untreated control for the three types of epithelia tested in this study.

Altogether, microscopy data evidenced an IL-13 dose-dependent increase of MUC5AC for all three types of epithelia, MucilAir™-Pool nasal, MucilAir™-Pool bronchial and SmallAir™-Pool using both analysis methods. Co-addition of Lebrikizumab\_ a commercially available IL-13 antagonist\_ to IL-13 efficiently prevented MUC5AC increase. Notably, MucilAir™-Pool Bronchial epithelium gives the best response to IL-13 exposure in term of MUC5AC production with a maximum value of 16 % MUC5AC area measured by IHC and 37255  $\mu\text{m}^3$  measured in 3D IF. Considering the fold increase of MUC5AC % area and volume, the bronchial epithelium was behind the bronchiolar tissue (Figure 7).

MUC5AC gene expression analysis also indicated a dose-dependent increase although the effect was moderate for MucilAir™-Pool Nasal and SmallAir™-Pool. These data are in line with microscopy results as bronchial epithelium displays the highest MUC5AC expression levels compared to the other tissue types. Based on all these data, we concluded that MucilAir™-Pool Bronchial and SmallAir™-Pool were the tissues of choice to study metaplasia. MucilAir™-Pool nasal also responded to IL-13 but the effect was moderate (Figure 6 and 7).

### Conclusions

- The new approach developed by Epithelix to quantify MUC5AC, based on combined immunofluorescence and 3D confocal microscopy picture analysis (Figure 6B) is operational as it led to results that were aligned with immunohistochemistry traditional method (See Figure 6A).
- MucilAir™-Pool Bronchial epithelium is the tissue of choice for studying metaplasia as it gave the highest response in term of MUC5AC production. SmallAir™-Pool is also a valuable epithelium for studying metaplasia in terms of MUC5AC fold change.

