

# Novel fully primary human airway epithelium-alveolar macrophages in vitro co-cultures models to study host pathogen interactions

Bernadett Boda, Carole Bertinetti, Ophélie Verbeke, Gowsinth Gunasingam, Song Huang and Samuel Constant

Epithelix Sàrl, 18 chemin des Aulx, CH-1228 Plan-les-Ouates, Geneva, Switzerland

## Introduction

Being the body's first line of defense against airborne pathogens like bacteria and viruses, the respiratory epithelium acts as a physical barrier as well as an effective mucociliary escalator. Moreover, the airway epithelium is also a powerful immunoregulator that orchestrates innate and adaptive immune responses during bacterial or viral infections. Here we established a new co-culture model using human pulmonary macrophages (M $\Phi$ ) and well characterized, standardized human airway epithelium in air-liquid interface culture, such as MucilAir<sup>TM</sup>, SmallAir<sup>TM</sup>.

# 1. Methods

1. 1 Isolation of pulmonary macrophages

1. 2 Characterization of pulmonary macrophages



Figure 1. Human pulmonary macrophages were isolated from lung provided by biobanks with informed consent and with approval from local ethics commissions.



Figure 2. Characterization of pulmonary macrophages. A. Cellular morphology of macrophages from two different donors (donor 1 and 2) using cytospin and Giemsa staining. B. Flow cytometry of pulmonary macrophages. The panel of macrophage markers was designed according to Bharat et al., 2016. Example from a healthy donor. (donor 2)



# 2. Results

2.1 Co-culture of airway epithelium and pulmonary macrophages







2.2. Phagocytic activity of the macrophages in co-culture







Figure 5. Nasal MucilAir<sup>TM</sup>-Pool-M $\Phi$  (donor 2) was infected with 10<sup>6</sup> CFU of MSSA (ATCC 29213) in 10  $\mu$ L. A. General schema for endpoint measurements. Co-culture was compared to infected MucilAir<sup>TM</sup>-Pool and infected MucilAir<sup>™</sup>-Pool treated with vehicle or appropriate antibiotic. Endpoints were assessed at 24 hours post infection (4 independent experiments, n=12, mean+SEM). B. Transepithelial electrical resistance (TEER). C. Lactate dehydrogenase (LDH) assay (Triton X-100: 100 % toxicity. D. Enumeration of bacteria at the apical side. One-way ANOVA with Dunnett's posttest, p<0.05 (GraphPad, Prism).

#### Prism).

### Summary

- > Pulmonary macrophages are CD45+/HLADR+ and show donor specific expression of CD206, CD169, CD11b, CD68 and CD14 markers
- Macrophages co-cultured with MucilAir<sup>TM</sup>/SmallAir<sup>TM</sup> are viable and have phagocytic activity in the co-culture
- > Macrophages in co-culture decrease the growth of Staphylococcus *aureus* (MSSA) by more than 1 log and prevent bacterium-induced toxicity
- > Macrophages in co-culture decrease the growth of Streptococcus pneumoniae (Sp19F) by more than 3 log and prevent bacteriuminduced increase of mucin secretion

#### Conclusion

Novel fully primary human airway epithelium-alveolar macrophages in vitro co-cultures models are promising research tools to study infectious diseases.