

## In vitro human airway epithelial platform for the development of novel anti-bacterial drugs

Laureen Jaupart, Christine Caul-Futy, Mendy Bouveret, Ophelie Verbeke, Carole Bertinetti, Mireille Caul-Futy, Samuel Constant Epithelix Sàrl, 18 chemin des Aulx, 1228 plan-les-Ouates, Geneva, Switzerland

Respiratory bacterial infections cause frequently mild to severe diseases worldwide. To develop new anti-bacterial drugs more predictive research tools are needed. We report herein the use of 3D epithelia made of primary human airway epithelial cells, MucilAir<sup>TM</sup>, for anti-bacterial drug screening. As proof-of-concept, typical disease-causing bacteria is used to infect human nasal epithelia reconstituted from a pool of 14 donors.

**Tested bacteria** 

1000

Pseudomonas aeruginosa (Pa), Staphylococcus aureus (Sa), Acinetobacter baumannii (Ab), Streptococcus pneumoniae (Sp), Haemophilus influenzae (Hi)

Inhibitor of bacterial growth used Meropenem (Mero), Amoxicillin (Amox) Cytomix for inflammation, Triton X-100 for cytotoxicity **Positive controls** Vehicle (VH) **Negative controls** 

 $\geq$ 

treated with alveolar macrophages.

B. Mucin secretion 24 hours post-

infection. Vehicle control was

compared to infected condition and

infected condition treated with

(mean±SEM). Student's t test

(Prism 6.0 GraphPad, \*p<0.05,

\*\*p<0.01, \*\*\*p<0.001, \*\*\*\*p<0.0001,

macrophages.

alveolar

ns: not significant).



appropriate antibiotic (mean±SEM).

Student's t test (Prism 6.0 GraphPad,

\*\*\*p<0.001,


50.000 primary

alveolar macrophages

\*\*p<0.01,

\*\*\*\*p<0.0001, ns: not significant).

\*p<0.05,



## CONCLUSION

Fingerprint of bacterial-specific effects on standardized in vitro nasal epithelium (MucilAir<sup>TM</sup> Pool-Nasal) are herein reported. Antibiotics like Meropenem inhibits the bacterial growth and abrogates its side effects. Similarly, macrophages in co-culture decrease the growth of Sp and prevent the bacterium-induced increase of mucin secretion. Additionnal end-point can be used such as mucociliary clearance (MCC), intratissular bacterial localisation, etc... These results suggest that MucilAir<sup>TM</sup> is a reliable tool for anti-bacterial drug development.